



# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g55630 (At1g55630)

<b>Product Code</b>	CSB-MP773700DOA
<b>Abbreviation</b>	At1g55630
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7X6A5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNSVIHYSTS VAVRKASRFL FTSRKFCNGS IGGDVTDNGT EEPLKITWES SEMDCEFDQE ENGEKISVRK RFMESTKLSA SRVLDLQGD CPGFNTKSAL DELNVSISGL LVREVLVGIL RTLSFDNKTR CAKLAYKFFV WCGGQENFRH TANCYHLLMK IFAECGEYKA MCRLIDEMIK DGYPTTACTF NLLICTCGEA GLARDVVEQF IKSCTFNYP YKHSYNAILH SLLGVKQYKL IDWVYEQMLE DGFTPDLTY NIVMFANFRL GKTDRLYRLL DEMVKDGFSP DLYTYNILLH HLATGNKPLA ALNLLNHMRE VGVEPGVIHF TTLIDGLSRA GKLEACKYFM DETVKVGCTP DVVCYTVMIT GYISGGELEK AEEMFKEMTE KGQLPNVFTY NSMIRGFCMA GKFKEACALL KEMESRGCNP NFVYVSTLVN NLKNAGKVLE AHEVVKDMVE KGHYVHLISK LKKYRRS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	At1g55630
<b>Protein Names</b>	Recommended name: Pentatricopeptide repeat-containing protein At1g55630
<b>Expression Region</b>	1-477
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.