



Recombinant *Oryza sativa* subsp. japonica Probable E3 ubiquitin-protein ligase BAH1-like 2 (Os07g0673200, LOC_Os07g47590)

Product Code	CSB-EP773711OFG
Abbreviation	Os07g0673200, LOC_Os07g47590
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XI73
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MKFAKKEYEY MKGMDEELPG VGLKRLKLL KRCRSDLQSH ENDGSSAGRC PGHCSVCDGS FFPSELLNEMS AVIGCFNEKA KKLLEHLAS GFKKYTMWFT SKGHKSHGAL IQQGKDLVTY AINAVAMRK ILKKYDKIHY SKQQQEFKAQ AQLHIEILQ SPWLCELMF YMNLRRSKKN NGAMELFGDC SLVFDDDKPT ISCNLFDSMR VDISLTCISIC LDTVFDPAVAL SCGHIYCYLC SCSAASVTIV DGLKSAERKS KCPLCRQAGV FPNVHLEDEL NMLLSYSCPE YWEKRIQMER VERVRLAKEH WESQCRAFLG M
Source	<i>E. coli</i>
Target Names	Os07g0673200
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase BAH1-like 2 EC=6.3.2.-
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.