



Recombinant *Oryza sativa* subsp. japonica Probable L-ascorbate peroxidase 7, chloroplastic (APX7)

Product Code	CSB-BP773712OFG
Abbreviation	APX7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XJ02
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	CMASASASA ASAAVASGAA ELKAAREDIR ELLKTTHCHP ILVRLGWHDS GTYDKNIKEW PQRGGANGSL RFDVELKHGA NAGLVNALKL VQPIKDKYPN ISYADLFQLA SATAIEEAGG PKIPMTYGRI DVTGPEQCPP EGKLPDAGPS APADHLRKVF YRMGLDDKEI VVLSGAHTLG RSRPERSGWG KPETKYTKNG PGAPGGQSWT AEWLKFDNSY FKEIKEKRDQ DLLVLPTDAA LFEDPTFKVY AEKYAEDQEA FFKDYAGAHA KLSNLGAKFN PPEGFTLDG
Source	Baculovirus
Target Names	APX7
Protein Names	Recommended name: Probable L-ascorbate peroxidase 7, chloroplastic EC=1.11.1.11 Alternative name(s): OsAPx07
Expression Region	72-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.