



Recombinant *Oryza sativa* subsp. japonica Cysteine protease 1 (CP1)

Product Code	CSB-EP773722OFG
Abbreviation	CP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XR52
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	LPDSVD WRDKGAVVAP VKNQGQCGSC WAFSAVAAVE GINKIVTGEL VSLSEQELVE CARNGQNSGC NGGIMDDAFA FIARNGGLDT EEDYPYTAMD GKCNLAKRSR KVVSIDGFED VPENDELSLQ KAVAHQPVSV AIDAGGREFQ LYDSGVFTGR CGTNLDHGVV AVGYGTDAAT GAAYWTVRNS WGPDWGENGY IRMERNVTAR TGKCGIAMMA SYPIKKG P
Source	E.coli
Target Names	CP1
Protein Names	Recommended name: Cysteine protease 1 EC= 3.4.22.- Alternative name(s): OsCP1
Expression Region	155-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.