



# Recombinant *Oryza sativa* subsp. *japonica* Serine/threonine-protein kinase SAPK8 (SAPK8)

<b>Product Code</b>	CSB-EP773735OFG-B
<b>Abbreviation</b>	SAPK8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7Y0B9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAAGAGAGA PDRAALTVGP GMDMPIMHDS DRYELVRDIG SGNFGVARLM RDRRTMELVA VKYIERGEKI DDNVQREIIN HRSLKHPNII RFKEVILTPT HLAIVMEYAS GGELFERICK NVRFSEDEAR YFFQQLISGV SYCHSMQVCH RDLKLENTLL DGSPAPRLKI CDFGYSKSSV LHSQPKSTVG TPAYIAPEVL LKKEYDGKTA DVWSCGVTLY VMVVGAYPFE DPEEPKNFRK TIQRILNVQY SIPENVDISP ECRHLISRIF VGDPSLRITI PEIRSHGWFL KNLPADLMDD DSMSSQYEEP DQPMQTMQI MQILTEATIP PACSRINHIL TDGLDLDDDM DDLSDSDID VDSSGEIVYA M
<b>Source</b>	E.coli
<b>Target Names</b>	SAPK8
<b>Protein Names</b>	Recommended name: Serine/threonine-protein kinase SAPK8 EC= 2.7.11.1 Alternative name(s): Osmotic stress/abscisic acid-activated protein kinase 8
<b>Expression Region</b>	1-371
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.