



Recombinant Human Bcl-2/adenovirus E1B 19 kDa-interacting protein 2-like protein (BNIPL)

Product Code	CSB-EP773788HU
Abbreviation	BNIPL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z465
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGTIQEAGKK TDVGVREIAE APELGAALRH GELELKEEWQ DEEFPRLLPE EAGTSEDPED PKGDSQAAAG TPSTLALCGQ RPMRKRLSAP ELRLSLTKGP GNDGASPTQS APSSPDGSSD LEIDELETPS DSEQLDGHE FEWEDELPA EGLGTSETAE RLGRGCMWDV TGEDGHHWRV FRMGPREQRV DMTVIEPYKK VLSHGGYHGD GLNAVILFAS CYLPRSSIPN YTYVMEHLFR YMVGTLELLV AENYLLVHLS GGTSRAQVPP LSWIRQCYRT LDRRLRKNLR ALVVVHATWY VKAFLALLRP FISSKFTRKI RFLDSL GELA QLISLDQVHI PEAVRQLDRD LHGSGGT
Source	E.coli
Target Names	BNIPL
Protein Names	Recommended name: Bcl-2/adenovirus E1B 19 kDa-interacting protein 2-like protein
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.