



Recombinant Human Amyloid beta A4 precursor protein-binding family B member 1-interacting protein (APBB1IP)

Product Code	CSB-BP773797HU
Abbreviation	APBB1IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z5R6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGESSEDIDQ MFSTLLGEMD LLTQSLGVDT LPPDPNPPR AEFNYSVGFK DLNESLNALE DQDLDALMAD LVADISEAEQ RTIQAQKESL QNQHHSASLQ ASIFSGAASL GYGTVNAATG ISQYEDDLPP PPADPVLDP LPPPPPEPLS QEEEEAQAKA DKIKLALAKL KEAKVKKLVV KVHMNDNSTK SLMVDERQLA RDVLDNLFK THCDCNVDWC LYEIYPELQI ERFFEDHENV VEVLSDWTRD TENKILFLEK EEKYAVFKNP QNFYLDNRGK KESKETNEKM NAKNKESLLE ESFCGTSIIV PELEGALYLK EDGKKS WKRR YFLLRASGIY YVPGKGTKTS RDLACFIQFE NVNIYYGTQH KMKYKAPTDY CFVLKHPQIQ KESQYIKYLC CDDTRTLNQW VMGIRIAKYG KTLYDNYQRA VAKAGLASRW TNLGTVNAAA PAQPSTGPKT GTTQPNGQIP QATHSVSAVL QEAQRHAETS KDKKPALGNH HDAVPRAPH APKSSLPPPP PVRSSDTSG SPATPLKAKG TGGGGLPAPP DDFLPPPPPP PPLDDPELPP PPPDFMEPPP DFVPPPPPSY AGIAGSELPP PPPPPPAPAP APVPDSARPP PAVAKRPPVP PKRQENPGHP GGAGGGEQDF MSDLMKALQK KRG NVS
Source	Baculovirus
Target Names	APBB1IP
Protein Names	Recommended name: Amyloid beta A4 precursor protein-binding family B member 1-interacting protein Alternative name(s): APBB1-interacting protein 1 Proline-rich EVH1 ligand 1 Short name= PREL-1 Proline-rich protein 73 Rap1-GTP-in
Expression Region	1-666
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.