



Recombinant Human Transcription factor AP-2-delta (TFAP2D)

Product Code	CSB-MP773802HU
Abbreviation	TFAP2D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z6R9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSTTFPGLVH DAEIRHDGSN SYRLMQLGCL ESVANSTVAY SSSSPLTYST TGTEFASPYF STNHQYTPLH HQSFHYEFQH SHPAVTPDAY SLNSLHHSQQ YYQQIHHGEP TDFINLHNAR ALKSSCLDEQ RRELGLDAY RRHDLSLMSH GSQYGMHPDQ RLLPGPSLGL AAAGADDLQG SVEAQCGLVL NGQGGVIRRG GTCVVNPTDL FCSVPGRLSL LSSTSKYKVT IAEVKRRLSP PECLNASLLG GILRRAKSKN GGRCLREKLD RLGLNLPAGR RKAANVTLT SLVEGEALHL ARDFGYTCET EPPAKAVGEH LARQHMEQKE QTARKKMILA TKQICKEFQD LLSQDRSPLG SSRPTPILDLDIQRHLTHFS LITHGFGTPA ICAALSTFQT VLSEMLNYLE KHTTHKNGGA ADGQGHANS EKAPLRKTSE AAVKEGKTEK TD
Source	Mammalian cell
Target Names	TFAP2D
Protein Names	Recommended name: Transcription factor AP-2-delta Short name= AP2-delta Alternative name(s): Activating enhancer-binding protein 2-delta Transcription factor AP-2-beta-like 1
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.