



Recombinant Human Patatin-like phospholipase domain-containing protein 5 (PNPLA5)

Product Code	CSB-MP773803HU
Abbreviation	PNPLA5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z6Z6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGFLEEEGRW NLSFSGAGYL GAHHVGATEC LRQRAPRLIQ GARRIYGSSS GALNAVSIVC GKSVDFFCCSH LLGMVGQLER LLSILHPAY APIEHVKQQL QDALPPDAHV LASQRLGISL TRWPDGRNFL VTDFATCDEL IQALVCTLYF PFYCGLIPPE FRGERYIDGA LSNLFPADC PSTITVSPFH GTVDICPQST SPNLHELNVF NFSFQISTEN FFLGLICLIP PSLEVVADNC RQGYLDALRF LERRGLTKEP VLWTLVSKEP PAPADGNWDA GCDQRWKGGGL SLNWKVPHVQ VKDVPNFEQL SPELEAALKK ACTRDPSRWA RFWHSRGPQV LTYLLLPTL PFEYIYFRSR RLVVWLPDVP ADLWWMQGLL RNMALEVFSR TKAQLLGPI SPPATRVLETS PLQPQIAPHR EELGPTHQA
Source	Mammalian cell
Target Names	PNPLA5
Protein Names	Recommended name: Patatin-like phospholipase domain-containing protein 5 EC= 3.1.1.- Alternative name(s): GS2-like protein
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.