



# Recombinant *Xenopus laevis* WD repeat domain phosphoinositide-interacting protein 2 (wipi2)

<b>Product Code</b>	CSB-YP773837XBE
<b>Abbreviation</b>	wipi2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7ZWU5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus laevis</i> (African clawed frog)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNLASQSVDA GAGQLLFANF NQDNTSLAVG SKSGYKFFSL SSVDKLEQIY ECTDTEVCI VERLFSSSLV AIVSLKAPRK LKVCHFCKGT EICNYSYSNT TLAVKLNRRQ LIVCLEESLY IHNIRDMKVL HTIRETPPNP SGLCSLSING ENCYLAYPGS ASIGEVQVFD TVNLRAANMI PAHDSPLAAL AFDASGTKLA TASEKGTVIR VFSIPEGQKL FEFRRGVKRC VSICSLAFSM DSIFLSASSN TETVHIFKLE TIKEPPEEP TSWTGYFGRV IMASTSYLPS QVTEMFNQGR AFATVRLPFC GHKNICALAT IQKISRLLVG AADGYLYIYN FDPQEGGECT LMKQHKLDGS MEPSSEILES SSHDRQVGAQ TYSATVTKTY PPPSPTRHAY ADDLGAVGGA SEEDEMGNLR LDEDNENPPM ILQTE
<b>Source</b>	Yeast
<b>Target Names</b>	wipi2
<b>Protein Names</b>	Recommended name: WD repeat domain phosphoinositide-interacting protein 2 Short name= WIPI-2
<b>Expression Region</b>	1-435
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.