



Recombinant Danio rerio Serine/threonine-protein phosphatase 2A regulatory subunit B' subunit gamma (ppp2r3c)

Product Code	CSB-BP773866DIL
Abbreviation	ppp2r3c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q803V3
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MANDITAHWK DLLRKRLASL KPDGRTEEEK KAESELSFSK YYTEWKGGEK GEDDSFKHIP RFYYRLPAED EVLMQKLREE SRAVFLQRKS RELLDNEELQ NLWFLLDKHQ VPPTTGDEAM ISYESFLKVG EKAGTKCKLF FTARVYAKLL HNDPYGRISI MQFFNYVMRK VWLHQTRIGL SLYDVAGQGY LRESLENYI LELIPTLPQL DGLEKSFYSF YVCTAVRKFF FFLDPLHTGK IKIQDILACS FLDDLLELRD EELSKESQES NWFSAPSALR VYGQYLNLDK DHNGMLSKEE LSRYGTGTLT SVFLDRVYQA CLTYDGEMDY KTYLDFVLAL ENRKEPAALQ YIFKLLDMEN KGYLNVFALN YFFRAIQEQM KIHGQEPVSF QDVKDEIFDM VKPKDPYKIT LQDLVNSGQG DTVSSILIDL NGFWTYENRE VLVANDTDSN AADLDDT
Source	Baculovirus
Target Names	ppp2r3c
Protein Names	Recommended name: Serine/threonine-protein phosphatase 2A regulatory subunit B' subunit gamma
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.