



# Recombinant Danio rerio Serine/threonine-protein phosphatase 2A regulatory subunit B' subunit gamma (ppp2r3c)

<b>Product Code</b>	CSB-MP773866DIL
<b>Abbreviation</b>	ppp2r3c
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q803V3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MANDITAHWK DLLRKRLASL KPDGRTEEEK KAESELSFSK YYTEWKGGEK GEDDSFKHIP RFYYRLPAED EVLMQKLREE SRAVFLQRKS RELLDNEELQ NLWFLLDKHQ VPPTTGDEAM ISYESFLKVG EKAGTKCKLF FTARVYAKLL HNDPYGRISI MQFFNYVMRK VWLHQTRIGL SLYDVAGQGY LRESLENYI LELIPTLPQL DGLEKSFYSF YVCTAVRKFF FFLDPLHTGK IKIQDILACS FLDDLLELRD EELSKESQES NWFSAPSALR VYGQYLNLDK DHNGMLSKEE LSRYGTGTLT SVFLDRVYQA CLTYDGEMDY KTYLDFVLAL ENRKEPAALQ YIFKLLDMEN KGYLNVFALN YFFRAIQEQM KIHGQEPVSF QDVKDEIFDM VKPKDPYKIT LQDLVNSGQG DTVSSILIDL NGFWTYENRE VLVANDTDSN AADLDDT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ppp2r3c
<b>Protein Names</b>	Recommended name: Serine/threonine-protein phosphatase 2A regulatory subunit B' subunit gamma
<b>Expression Region</b>	1-457
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.