



Recombinant *Xenopus laevis* Zinc finger protein 36, C3H1 type-like 2-B (zfp36l2-B)

Product Code	CSB-MP773869XBE
Abbreviation	zfp36l2-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q805B4
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MSTLLSAFY DIDLLYKNEK ALNNLALSTM LDKKAVGSPV SSTNSNLFPG FLRRHSATNL QALSGSTNLA KFCPNNNNNP LKDPAVSSTA LLNRENKFRD RSFSENGERS QHLLHLQQQQ QQQKAGAQVN STRYKTELCR PFEENGACKY GEKCQFAHGF HELRSLTRHP KYKTELCRTF HTIGFCPYGP RCHFIHNAEE RRQAPGAGER PKLHHSLSFS GFPNHSLDSP LLESPTSRT PPQSSGSLYC QELLQLNNNN PCANNAFTFS GQELGLIAPL AIHTQNQSYC RQPCSSPPLS FQPLRRVSES PVFDAPPSPD DSLSDRDSYL SGSLSSGSL GSDSPTLDSN RRLPIFSRLS ISDD
Source	Mammalian cell
Target Names	zfp36l2-B
Protein Names	Recommended name: Zinc finger protein 36, C3H1 type-like 2-B Alternative name(s): CCCH zinc finger protein 3-B Short name= XC3H-3b
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.