



Recombinant Oligoribonuclease (orn)

| | |
|--------------------------|---|
| Product Code | CSB-YP774476SZB |
| Abbreviation | orn |
| Storage | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |
| Uniprot No. | Q83IJ9 |
| Product Type | Recombinant Protein |
| Immunogen Species | Shigella flexneri |
| Purity | >85% (SDS-PAGE) |
| Sequence | SANENNLIW IDLEMTGLDP ERDRIIEIAT LVTDANLNIL AEGPTIAVHQ SDEQLALMDD WNVRTHTASG LVERVKASTM GDREAELATL EFLKQWVPAG KSPICGNSIG QDRRFLFKYM PELEVYFHYR YLDVSTLKEK ARRWKPEILD GFTKQGTHQA MDDIRESVAE LAYYREHFIK L |
| Source | Yeast |
| Target Names | orn |
| Protein Names | Recommended name: Oligoribonuclease EC= 3.1.-.- |
| Expression Region | 2-181 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |