



Recombinant Putative multidrug resistance protein MdtA (mdtA)

Product Code	CSB-BP774494SZB
Abbreviation	mdtA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83KI5
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MLDNLIGARY LTGLGTITAA NVTVRSRVD GQLMALHFQE GQQVKAGDLL AEIDPSQFKV ALAQTQGQLA KDKATLANAR RDLARYQQLA KTNLVSQRQEL DAQQALVSET EGTIKADEAS VASAQLQLDW SRITAPVDGR VGLKQVDVGN QISSGDTTGI VVITQTHPID LVFTLPESDI ATVVQAQKAG KPLMVEAWDR TNSKKLSEGT LLSLDNQIDA TTGTIKVKAR FNNQDDALFP NQFVNARMLV DTEQNAVVIP TAALQMGNEG HFVWVLNSEN KVSKHLVTPG IQDSQKVVIR AGISAGDRVV TDGIDRLTEG AKVEVVEAQS ATTPEEKATS REYAKKGARS
Source	Baculovirus
Target Names	mdtA
Protein Names	Recommended name: Putative multidrug resistance protein MdtA Alternative name(s): Multidrug transporter MdtA
Expression Region	1-350
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.