



Recombinant Deferrochelataase/peroxidase EfeB (efeB)

Product Code	CSB-BP774503SZB
Abbreviation	efeB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83LK4
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	QKTQSAP GTLSPDARNE KQPFYGEHQA GILTPQQAAM MLVAFDVLAS DKADLERLFR LLTQRFAFLS QGGAAPETPN PRLPPLDSGI LGGYIAPDNL TITLSVGHS� FDERFGLAPQ MPKKLQKMTR FPNDSLDAAL CHGDVLLQIC ANTQDTVNHA LRDIIKHTPD LLSVRWKREG FISHAARSK GKETPINLLG FKDGTANPDS QNDKLMQKVW WVTADQQEPA WTIGGSYQAV RLIQFRVEFW DRTPLKEQQT IFGRDKQTGA PLGMQHEHDV PDYASDPEGK VIALDSHIRL ANPRTAESES SLMLRRGYSY SLGVTNSGQL DMGLLFVCYQ HDLEKGFLTIV QKRLNGEAL EYVKPIGGGY FFSLPGVKDA NDYLGRALLQ V
Source	Baculovirus
Target Names	efeB
Protein Names	Recommended name: Deferrochelataase/peroxidase EfeB EC= 1.11.1.-
Expression Region	34-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.