



Recombinant Taurine import ATP-binding protein TauB (tauB)

Product Code	CSB-MP774508SZB
Abbreviation	tauB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83MA0
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MLQISHLYAD YGGKPALEDI NLTLESGELL VVLGPSGCGK TTLNLIAGF VPYQHGSIQ L AGKRIEGPGA ERGVVFQNEG LLPWRNVQDN VAFGLQLAGI EKMQRLEIAH QMLKKVGLEG AEKRYIWQLS DGQRQRVGIA RALANPQLL LLDEPGALD AFTRDQMOTL LLKLWQETGK QVLLITHDIE EAVFMATELV LLSSGPGRVL ERLRLNFARR FVAGESSRSI KSDPQFIAMR EYVLSRVFEQ REAFS
Source	Mammalian cell
Target Names	tauB
Protein Names	Recommended name: Taurine import ATP-binding protein TauB EC= 3.6.3.36
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.