



Recombinant 4-hydroxy-2-oxo-heptane-1,7-dioate aldolase (hpcH)

Product Code	CSB-YP774524SZB
Abbreviation	hpcH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83P19
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MENSFKAALK AGRPQIGLWL GLSSSYS AEL LAGPGFDWLL IDGEHAPNNV QTVLTQLQAI APYPSQPVVR PSWNDPVQIK QLLDVGTQTL LVPMVQNADE ARKAVRATRY PPAGIRGVGS ALARASRWRN IPDY LQKAND QMCVLVQIET REAMKNLPQI LDVEGV DGVF IGPADLSADM GYAGNPQHPE VQAAIEQAIV QIREAGKAPG ILIANEQLAK RYLELGALFV AVGVDTTLLA RAAEALAARF GNSSSISMAK QNNNSVY
Source	Yeast
Target Names	hpcH
Protein Names	Recommended name: 4-hydroxy-2-oxo-heptane-1,7-dioate aldolase EC= 4.1.2.n4 Alternative name(s): 2,4-dihydroxyhept-2-ene-1,7-dioic acid aldolase Short name= HHED aldolase 4-hydroxy-2-ketoheptane-1,7-dioate aldolase Short name=
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.