



# Recombinant Arabidopsis thaliana Putative DNA repair protein RAD23-1 (RAD23-1)

<b>Product Code</b>	CSB-YP774596DOA
<b>Abbreviation</b>	RAD23-1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q84L33
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKLTVKTLKG SHFEIRVLPS DTIMAVKKN I EDSQGKDNYP CGQQLLIHNG KVLKDETSLV ENKVTEEGFL VVMLSKSKSG GSAGQASVQT SSVSQPVSAT TSSTKPAAPS TTQSSPVPAS PIPAQEQPAA QTDTYGQAAS TLVSGSSLEQ MVQQIMEMGG GSWDKETVTR ALRAAYNNPE RAVDYLYSGI PQTAEVAVPV PEAQIAGSGA APVAPASGGP NSSPLDLFPQ ETVAAAGSGD LGTLEFLRNN DQFQQLRTMV HSNPQILQPM LQELGKQNPQ LLRLIQENQA EFLQLVNEPY EGSDGEGDMF DQPEQEMPHA INVTPAEQEA IQRLEAMGFD RALVIEAFLA CDRNEELAN YLLENSGDFE D
<b>Source</b>	Yeast
<b>Target Names</b>	RAD23B
<b>Protein Names</b>	Recommended name: Putative DNA repair protein RAD23-1 Alternative name(s): RAD23-like protein 1 Short name= AtRAD23-1
<b>Expression Region</b>	1-371
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.