



# Recombinant *Oryza sativa* subsp. japonica Cyclin-dependent kinase F-3 (CDKF-3)

<b>Product Code</b>	CSB-BP774609OFG
<b>Abbreviation</b>	CDKF-3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q84SN3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MERYKVIREI GDGTTCGNVFR AYNTEETNEIV AVKKMKRKFF QWEECISLRE VKALQKLNHP NIVKLKEVTM ENHELFFIFE NMECNLYDVI RERQAAFSEE EIRNFMVQIL QGLAYMHNNG YFHRDLKPEN LLVTDGTVKI ADFGLAREVS SSPPYTDYVS TRWYRAPEVL LQSSAYTPAI DMWAVGAILA ELFTLSPLFP GGSETDQLYK ICAVLGTPDH TVWPEGMNL P R S S S F N F F Q I P P R N L W E L I P NATLEAIDLI QQLCSWDPRR RPTAEQSLQH PFFNVGNWVP RPLHASHTKT IETRPNPRLE LNLWDFGTEP EDNYLDLTL LKPSFPGTDF SNNVPEHTKE EILLYPGFEN PPVQSGFWPL VASDRPMGDV PAMSSWPQAY VVDGQATLPA VGFSGSPFGL SPLQPNLFEN RSFATPIRQV NFF
<b>Source</b>	Baculovirus
<b>Target Names</b>	CDKF-3
<b>Protein Names</b>	Recommended name: Cyclin-dependent kinase F-3 Short name= CDKF;3 EC= 2.7.11.22 EC= 2.7.11.23 Alternative name(s): Serine/threonine-protein kinase MHK-like protein 1
<b>Expression Region</b>	1-433
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.