



# Recombinant Human MOB kinase activator 3B (MOB3B)

<b>Product Code</b>	CSB-YP774795HU
<b>Abbreviation</b>	MOB3B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86TA1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSIALKQVFN KDKTFRPKRK FEPGTQRFEL HKRAQASLNS GVDLKAHVQL PSGEDQNDWV AVHVVDFFNR INLIYGTICE FCTERTCPVM SGGPKYEYRW QDDLKYYKPT ALPAPQYMN LMDWIEVQIN NEEIFPTCVG VPFKPNFLQI CKKILCRLFR VVHVYIHHF DRVIVMGAEA HVNTCYKHFY YFVTEMNLID RKELEPLKEM TSRMCH
<b>Source</b>	Yeast
<b>Target Names</b>	MOB3B
<b>Protein Names</b>	Recommended name: MOB kinase activator 3B Alternative name(s): Mob1 homolog 2b Mps one binder kinase activator-like 2B
<b>Expression Region</b>	1-216
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein shares similarity with the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. This gene is located on the opposite strand as the interferon kappa precursor (IFNK) gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.