



Recombinant Human Protein FAM131B (FAM131B)

Product Code	CSB-MP774828HU
Abbreviation	FAM131B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86XD5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDSTSSLHGS SLHRPSTEQT RTDFSWDGIN LSMEDTTSIL PKLKRNSNAY GIGALAKSSF SGISRSMKDH VTKPTAMGQG RVAHMIEWQG WGKTPAVQPQ HSHESVRRDT DAYSDLSDGE KEARFLAGVM EQFAISEATL MAWSSMDGED MSVNSTQEPL GCNYSYDNYQE LMDSQDALAQ APMDGWPHSY VSQGMVCLGS SDAWEASDQS LIASPATGSY LGPAFDDSQP SLHEMGPSQP ASGYSALEPP PLLGGDTDWA PGVGAVDLAR GPAEEEEKRPL APEEEEDAGC RDLESLSPRE DPEMSTALSR KVSDVTSSGV QSFDEEEGEA NN
Source	Mammalian cell
Target Names	FAM131B
Protein Names	Recommended name: Protein FAM131B
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.