



Recombinant Human TIR domain-containing adapter molecule 2 (TICAM2)

Product Code	CSB-BP774830HU
Abbreviation	TICAM2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86XR7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GIGKSKINS CPLSLSWGKR HSVDTSPGYH EDSKKS EDDL SLCNVAEHSN TTEGPTGKQE GAQSVEEMFE EEAE EEFVLK FVILHAEDDT DEALRVQNL LQDDFGIKPGI IFAEMPCGRQ HLQNLDDAVN GSAWTILLT ENFLRDTWCN FQFYTSLMNS VNRQHKYNSV IPMRPLNNPL PRERTPFALQ TINALEEESR GFPTQVERIF QESVYKTQQT IWKETRNMVQ RQFIA
Source	Baculovirus
Target Names	TICAM2
Protein Names	Recommended name: TIR domain-containing adapter molecule 2 Short name= TICAM-2 Alternative name(s): Putative NF-kappa-B-activating protein 502 TRIF-related adapter molecule Toll-like receptor adaptor protein 3 Toll/interleukin-1
Expression Region	2-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.