



Recombinant tRNA dimethylallyltransferase (miaA)

Product Code	CSB-BP774954VFE
Abbreviation	miaA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q87L06
Product Type	Recombinant Protein
Immunogen Species	Vibrio parahaemolyticus serotype O3:K6 (strain RIMD 2210633)
Purity	>85% (SDS-PAGE)
Sequence	MTDKLPLALF LMGPTASGKT ELAIRLRQKF PVEIISVDSA LIYKGM DIGT AKPDERELSL APHRLIDILD PSESYSAADF RRDALQAMDD IVAEGKIPLL VGGTMLYYKA LLEGLSPLPA ANPEIRQQIE QEALTKGWSV LHDELKEIDP VSAARIHPND PQRLSRALEV YRISGKTLTE LTQTKGESLP YRVKQFAIAP KDRAELHRRRI ELRFDKMMEA GFEEEMKALY ARKDLHPDLP SIRCVGYRQM WEYLDGECTR DEAVFRGVCA TRQLAKRQIT WLRSWDDLTW LDSDNIEQAL ETMSEAIASD
Source	Baculovirus
Target Names	miaA
Protein Names	Recommended name: tRNA dimethylallyltransferase EC= 2.5.1.75 Alternative name(s): Dimethylallyl diphosphate:tRNA dimethylallyltransferase Short name= DMAPP:tRNA dimethylallyltransferase Short name= DMATase Isopentenyl-diph
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.