



Recombinant *Oryza sativa* subsp. *japonica* CBL-interacting protein kinase 2 (CIPK2)

Product Code	CSB-MP800014OFG
Abbreviation	CIPK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7X996
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAEQRGNMLM KKYEMGKLLG QGTFKAVYHA RNTETSESVA IKMIDKEKVL KGGLMDQIKR EISVMKLV RH PNIVQLYEVM ATKTKIYFVL EHVKG GELFN KVQRGRLKED AARKYFQQLI CAVDFCHSRG VYHRDLKPEN LLLDENS NLK VSDFGLSALA DCKRQDGLLH TTCGTPAYVA PEVINRRGYD GAKADIWSCG VILFVLLAGY LPFHDKNLMD MYKKIGKAEF KCPSWFNTDV RRLLLRILD P NPSTRISMDK IMENPWFRKG L DAKLLRYNL QPKDAIPVDM STDFDSFNSA PTLEKKPSNL NAFDIISLST GLDLSGMFEE SDKKESKFTS TSTASTIISK IEDIAKGLRL KLTKKD GLL KMEGSKPGRK GVMGIDAEIF EVTPNFHLVE LKKTNGDTLE YRKVLNQEMR PALKDIVWAW QGEQPKQQQQ PTC
Source	Mammalian cell
Target Names	CIPK2
Protein Names	Recommended name: CBL-interacting protein kinase 2 EC= 2.7.11.1 Alternative name(s): OsCIPK02
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.