



Recombinant *Oryza sativa* subsp. *japonica* Myb-related protein Myb4 (MYB4)

Product Code	CSB-MP800018OFG
Abbreviation	MYB4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XBH4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MGRAPCCEKM GLKKGPTWPE EDKVLVAHIQ RHGHGNWRAL PKQAGLLRCG KSCRLRWINY LRPDIKRGNF SKEEEDTIIH LHELLGNRWS AIAARLPGRT DNEIKNVWHT HLKKRLDAPA QGGHVAASGG KKHKPKPSAK KPAAAAAAPP ASPERSASSS VTESSMASSV AEEHGNAGIS SASASVCAKE ESSFTSASEE FQIDDSFWSE TLSMPLDGYD VSMEPGDAFV APPSADDMDY WLGVMESGE AQDLPQI
Source	Mammalian cell
Target Names	MYB4
Protein Names	Recommended name: Myb-related protein Myb4 Alternative name(s): OsMyb4 Transcription factor RLTR1
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.