



Recombinant *Oryza sativa* subsp. *japonica* SPX domain-containing protein 3 (SPX3)

Product Code	CSB-EP800021OFG-B
Abbreviation	SPX3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XEY9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MKFGKRLKKQ VEESLPEWRD KFLAYKRLKK LVRLVSSSSG DVGGGGGGGEA EFVRLLDGEV DRINAFFLEQ EEEFVIRQRE LQETVEKVAG GGGGGRRPAA AEMRRVRKEI VDLHGEMVLL LNYSVNYTG LAKILKKYDK RTGRLLRLPF IEKVLQRPF TTELISRLVR DCEATMEAIF TSSVATTAMA GDRRTWKGCS GDAGMAPMAD QQGIFRNTVA ALATMKELRS GSSTYGRFSL PPMAAPASPE SDVLQSIRSD PHLKENGRIP SLQFFYA
Source	E.coli
Target Names	SPX3
Protein Names	Recommended name: SPX domain-containing protein 3 Alternative name(s): Protein SPX DOMAIN GENE 3 Short name= OsSPX3
Expression Region	1-277
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.