



Recombinant *Oryza sativa* subsp. japonica

Probable protein phosphatase 2C 38

(Os04g0321800, LOC_Os04g25570)

Product Code	CSB-YP800043OFG
Abbreviation	Os04g0321800, LOC_Os04g25570
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XW27
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVAVTGGRPP GLQDAPGAPP PAPAAEAVPS RPLARDATYG GRVYGGVGGG GCCLEFLDCV LRAMGVATPA EIMPPADFRW AARPMRRRRR GGSSSSSSSP RDREPRDGRI AANGASAAAS LYTMRGNKGV NQDAMLVWEN FCSKEDTIFC GVFDGHGPGY HLVSKRVRDL LPIKLSANLG RDGHKETSTN IVTSSMTEGG GTERMDRDTE TPLGTEENGD YPEMFAALRT SLLRAFVMD RDLKFHKTID SVFSGTTAVT VIKQGHDLII GNLGDSRAVL GTRDEYDQFF AVQLTVDLKP TIPSEARIR ERSGRIFSLP DEPDVARVWL PKYNMPGLAM ARAFADFCLK DYGLISMPDV SYHRITEKDE FVVLATDGVW DVLSNSEVVS IVSQAKSEAS AARFVVESAQ RAWRTRFPTS KIDDCAVVCL FLNTDARNKP PGSGIKDLAN AIELGGGNLS
Source	Yeast
Target Names	Os04g0321800
Protein Names	Recommended name: Probable protein phosphatase 2C 38 Short name= OsPP2C38 EC= 3.1.3.16
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.