



Recombinant *Oryza sativa* subsp. japonica Zinc finger A20 and AN1 domain-containing stress-associated protein 9 (SAP9)

Product Code	CSB-YP800049OFG
Abbreviation	SAP9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Y1W9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAQESWKNES EETVHTPEAP ILCVNNGGFF GSSMTNNMCS KCYRDFVKVT TMAAPVVEKK AFTPASSSKT PLEPAKPDEV PAAAVEDKQA AQEPPKPPSN RCLSCRKKVG LTGFQCRCGG TFCSTHRYTE AHDCTFDYKK AGRDQIAKQN PVVIAEKINK I
Source	Yeast
Target Names	SAP9
Protein Names	Recommended name: Zinc finger A20 and AN1 domain-containing stress-associated protein 9 Short name= OsSAP9
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.