



Recombinant Human PTB-containing, cubilin and LRP1-interacting protein (PID1)

Product Code	CSB-EP800096HU-B
Abbreviation	PID1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z2X4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MWQPATERLQVTYLGKVSTTGMQFLSGCTEKPVIELWKKHTLAREDVFPANA LLEIRPFQ VWLHHLDDHKGEATVHMDTFQVARIAYCTADHNVSPNIFAWVYREINDDLSYQ MDCHAVEC ESKLEAKKLAHAMMEAFRKTFFHSMKSDGRIHSNSSSSEEVSQELESDDG
Source	E.coli
Target Names	PID1
Protein Names	Recommended name: PTB-containing, cubilin and LRP1-interacting protein Short name= P-CLI1 Alternative name(s): Phosphotyrosine interaction domain-containing protein 1 Protein NYGGF4
Expression Region	1-168
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.