



Recombinant Human Vitelline membrane outer layer protein 1 homolog (VMO1)

Product Code	CSB-EP800218HU
Abbreviation	VMO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z5L0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	QTDGRN GYTAVIEVTS GGPWGDWAWP EMCPDGGFFAS GFSLKVEPPQ GIPGDDTALN GIRLHCARGN VLGNTHTVVEV QSGSWGGEWSE PLWCRGGAYL VAFSLRVEAP TTLGDNTAAN NVRFRCSDGE ELQGPGLSWG DFGDWSDHCP KGACGLQTKI QGPRGLGDDT ALNDARLFCC RS
Source	E.coli
Target Names	VMO1
Protein Names	Recommended name: Vitelline membrane outer layer protein 1 homolog
Expression Region	25-202
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.