



Recombinant Human Sprouty-related, EVH1 domain-containing protein 2 (SPRED2)

Product Code	CSB-BP800221HU
Abbreviation	SPRED2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z698
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTEETHPDDD SYIVRVKAVV MTRDDSSGGW FPQEGGGISR VGVCKVMHPE GNRSGLFIH GERQKDKLVV LECYVRKDLV YTKANPTFHH WKVDNRKFGL TFQSPADARA FDRGVRKAIE DLIEGSTTSS STIHNEAELG DDDVFATTATD SSSNSSQKRE QPTRTISSPT SCEHRRIYTL GHLHDSYPTD HYHLDQPMPR PYRQVSFPDD DEEIVRINPR EKIWMGTGYED YRHAPVRGKY PDPSEDADSS YVRFKGEVP KHDYNYPYVD SDFGLGEDP KGRGGSVIKT QPSRGKSRRR KEDGERSRCV YCRDMFNHEE NRRGHCQDAP DSVRTCIRRV SCMWCADSM L YHCMSDPEGD YTDPCSCDTS DEKFCLRWMA LIALSFLAPC MCCYLPLRAC YHCGVMCRCC GGKHKAAA
Source	Baculovirus
Target Names	SPRED2
Protein Names	Recommended name: Sprouty-related, EVH1 domain-containing protein 2 Short name= Spred-2
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.