



Recombinant Human YTH domain-containing family protein 3 (YTHDF3)

Product Code	CSB-MP800225HU
Abbreviation	YTHDF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z739
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SATSVDQRP KGQGNKVSQ NGSIHQKDAV NDDDFEPYLS SQTNQSNSYP PMSDPYMPSY YAPSIGFPYS LGAAWSTAG DQPMPYLTTY GQMSNGEHY IPDGVFSQPG ALGNTPPFLG QHGFNFFPGN ADFSTWGTSG SQGQSTQSSA YSSSYGYPPS SLGRAITDGQ AGFGNDTLSK VPGISSIEQG MTGLKIGGDL TAAVTKTVGT ALSSSGMTSI ATNSVPPVSS AAPKPTSWAA IARKPAKPQP KLKPKGNVGI GGSAPPPPI KHNMNIGTWD EKGSVVKAPP TQPVLPPQTI IQQPQLIQP PPLVQSQLPQ QQPQPQPQQ QQGPQQAQP HQVQPQQQL QNRWVAPRNR GAGFNQNGA GSENFGLGVV PVSASPSSVE VHPVLEKLKA INNYNPKDFD WNLKNGRVFI IKSYSEDDIH RSIKYSIWCS TEHGNKRLDA AYRSLNGKGP LYLFLSVNGS GHFCGVAEMK SVVDYNAYAG VWSQDKWKGK FEVKWIFVKD VPNNQLRHIR LENNDNKPVT NSRDTQEVPL EKAKQVLKII ATFKHTTSIF DDFAHYEKRQ EEEEAMRRER NRNKQ
Source	Mammalian cell
Target Names	YTHDF3
Protein Names	Recommended name: YTH domain family protein 3
Expression Region	2-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.