



Recombinant Danio rerio H/ACA ribonucleoprotein complex subunit 1 (gar1)

Product Code	CSB-MP800251DIL
Abbreviation	gar1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZVE0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSFRGGGGGR GGGFNRGGGG GRGGGFGGGR GGGFGGGRGG GFGGGRGGRG GFNRNQDYGP PEYVVALGEF MHPCEDEIVC KCVTEENKVP YFNAPVYLEN KEQIGKVDEI FGQLRDFYFS VKLSDNMKAS SFKKLQKFYI DPMKLLPLQR FLPRPPGEKG PPRGGRGGGG GRGGRGGGFR GGRGANGGGR GGFGGRRGGGF GGRGGGGGGF RGGRGGGGGR GFRGGR
Source	Mammalian cell
Target Names	gar1
Protein Names	Recommended name: H/ACA ribonucleoprotein complex subunit 1 Alternative name(s): Nucleolar protein family A member 1 snoRNP protein GAR1
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.