



Recombinant Danio rerio Cyclin-L1 (ccn1)

Product Code	CSB-YP800255DIL
Abbreviation	ccn1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZVX0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSLGMLSPHL NTPPPNNQGI LIGDKVYSEV FLAIDNSIIP EDRLSTTPSM LDGLDHETET DLRILGCERI QSAGILLRLP QVAMATGQVI FQRFFFSKSF VKHNFEIVAM ACVNLASKIE ESPRRVRDVI NVFHHLKQGK GKKSTPLILD QNYINTKNQV IKAERRILKE LGFCVHVKHP HKIIVMYLQV LECEKNQMLV QTAWNYMNDALRTSAFVRFE PETIACACIY LAARVLQIPL PSKPHWFLLF GATKEDIKEI CINTMKLYSR EKPHSEQLER QVEKRKIFLE EARLKARGQN PNGTPALASI NGFSPASKPS SPRDVKMDDK SPNSKLKEPE NRQLFAKSPL NGSIKKEDGK VFQNGKNHSR SRSRSTSRSP HRHRRSHSGT YSSHSSHSPS PRQKARRPSP ISQLRTDRDR PSETSRHSNK RRRSRRSRS NSRERVRDRD HIKHKQERSG SGHHWDHRDR ERDRSRDHGR NKRQSRSHSG HSHSRHRR
Source	Yeast
Target Names	ccn1
Protein Names	Recommended name: Cyclin-L1
Expression Region	1-498
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.