



Recombinant Danio rerio 60S ribosomal protein L18a (rpl18a)

Product Code	CSB-MP800258DIL
Abbreviation	rpl18a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZWJ4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MKASGTLREY KVVGRLLPSV KNPAPPLYRM RIFAPNHVVA KSRFWYFVSQ LRKMKKASGE IVYCGLVFEK SPLKVKNFGI WLYDSSRS GT HNMRYREYRDL TTSGAVTQCY RDMGARHRAR AHAIQIMKVQ VIAANKCRRP AIKQFHDSKI KFPLPHRVLR RQHKPRFTTK RPNTFF
Source	Mammalian cell
Target Names	rpl18a
Protein Names	Recommended name: 60S ribosomal protein L18a
Expression Region	1-176
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.