



# Recombinant Mouse Protein deltex-3 (Dtx3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP800337MO  |
| <b>Abbreviation</b>      | Dtx3  |
| <b>Storage</b>           | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>   |
| <b>Uniprot No.</b>       | Q80V91  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSFVLSRMAA CGGSCKNKVT VSKPVWDFLS KETPARLARL REEHRVSILI<br>DGETSDIYVL QLSPQGPPPA PPNGLYLARK ALKGLLKEAE KELKKAQRQG<br>ELMGCLALGG GGEHPELHRP GPPPLRAAPL LPPGARGLPP PPPPLPPPLP<br>PRLREDAEEQ ETTCPICLGE IQNAKTLEKC RHSFCEGCIT RALQVKKACP<br>MCGRFYGQLV GNQPQNGRML VSKDATLLLP SYEKYGTIVI QYVFPPGVQG<br>AEHPNPGVRY PGTTRVAYLP DCPEGNKVLT LFRKAFDQRL TFTIGTSMTT<br>GRPNVITWND IHHKTSCTGG PQLFGYPDPT YLTRVQEELR AKGITDD |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | Dtx3  |
| <b>Protein Names</b>     | Recommended name: Protein deltex-3 Short name= Deltex3 Short name=<br>mDTX3   |
| <b>Expression Region</b> | 1-347   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.     |
| <b>Shelf Life</b>        | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>   |