



Recombinant Mouse Phosphatidylinositol 5-phosphate 4-kinase type-2 beta (Pip4k2b)

| | |
|--------------------------|--|
| Product Code | CSB-EP800358MO |
| Abbreviation | Pip4k2b |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q80XI4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSSNCTSTTA VAVAPLSASK TKTKKKHFVC QKVKLFRASE PILSVLMWGV NHTINELSNV PVPVMLMPDD FKAYSKIKVD NHLFNKENLP SRFKFKEYCP MVFRNLRERF GIDDQDYQNS VTRSAPINSQ SQRGCGTRFL TTYDRRFVIK TVSSEDVAEM HNILKKYHQF IVECHGNTLL PQFLGMYRLT VDGVETYMVV TRNVFSHRLT VHRKYDLKGS TVAREASDKE KAKDLPTFKD NDFLNEGQKL HVGEEESKKNF LEKLRDVEF LAQLKIMDYS LLVGIHDVDR AEQEEMEVEE RAEEEECEND GVGGSLLCSY GTPPDSPGNL LSFPRFFGPG EFDPSVDVYA MKSHEAPKK EVYFMAIIDI LTPYDAKKKA AHAAKTVKHG AGAEISTVNP EQYSKRFNEF MSNILT |
| Source | E.coli |
| Target Names | Pip4k2b |
| Protein Names | Recommended name: Phosphatidylinositol 5-phosphate 4-kinase type-2 beta EC= 2.7.1.149 Alternative name(s): 1-phosphatidylinositol 5-phosphate 4-kinase 2-beta Diphosphoinositide kinase 2-beta Phosphatidylinositol 5-phosphate 4-kinase |
| Expression Region | 1-416 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Target Details | This protein catalyzes the phosphorylation of phosphatidylinositol-5-phosphate on the fourth hydroxyl of the myo-inositol ring to form phosphatidylinositol-5,4-bisphosphate. This gene is a member of the phosphatidylinositol-5-phosphate 4-kinase family. The encoded protein sequence does not show similarity to other kinases, but the protein does exhibit kinase activity. Additionally, the encoded protein interacts with p55 TNF receptor. |

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.