



Recombinant Acetyl esterase (aes)

Product Code	CSB-EP800935SZB
Abbreviation	aes
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q83M39
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	≥85% (SDS-PAGE)
Sequence	MKPENKLPVL DLISAEMKTV VNTLQPDLPS WPATGTIAEQ RQYYTLERRF WNAGAPEMAT RAYMVPTKYG QVETRLFCPQ PDSPATLFYL HGGGFILGNL DTHDRIMRLL ASYSQCTVIG INYTLSPEAR FPQAIEEIVA ACCYFHQQAE DYQINMSRIG FAGDSAGAML ALASALWLRD KQIDCGKIAG VLLWYGLYGL RDSVTRLLG GVWDGLTQQD LQMYEEAYLS NDADRESPYY CLFNNDLTRE VPPCFIAGAE FDPDLLDSRL LYQTLAAHQQ PCEFKLYPGT LHAFLHYSRM MKTADALRD GAQFFTAQL
Source	E.coli
Target Names	aes
Protein Names	Recommended name: Acetyl esterase EC= 3.1.1.-
Expression Region	1-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>