



Recombinant HTH-type transcriptional regulator frIR (frIR)

Product Code	CSB-MP800963SZB
Abbreviation	frIR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83PX4
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MSATDRYSHQ LLYATVRQRL LDDIAQGVYQ AGQQIPTENE LCTQYNVSRI TIRKAISDLV ADGVLIRWQG KGTFVQSQKV ENALLTVSGF TDFGVSQGKS TKEKVIEQER VSAAPFCEKL NIPGNSEVFH LCRVMYLDKE PLFIDSSWIP LSRYPDFDEI YVEGSSTYQL FQERFDTRVV SDKKTIDIFA ATRPQAKWLK CELGEPLFRI SKIAFDQNDK PVHVSELCR ANRITLTIDN KRH
Source	Mammalian cell
Target Names	frIR
Protein Names	Recommended name: HTH-type transcriptional regulator frIR
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.