



Recombinant Arabidopsis thaliana E3 ubiquitin-protein ligase SINA-like 10 (At5g37930)

Product Code	CSB-BP801017DOA
Abbreviation	At5g37930
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84K34
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MARFSVCGGD DGEGPSNNNH QSRKRQRLPS IDENEEDAET SDAGSSGEED EDETQNQGMR PESEDRGSTS DDSREVVIE ERRFGKFNVS QSSSSKDSP LSVTLDPDV LDCPICCEPL KIPIFQCDNG HLACTLCCTK VRNRCPSCTL PIGYVRCRAM EKVIEASRVS CLNAKYGCKE STSYGNRFSH EQVCVFTPCS CPILDCHYTG YYKDLNNHVR AEHKDDLISF VWNTRLTISL DLNKKTTILQ EENDGHVIVV QVFRALHAVY VSVSCIAPLT PGVGRLSCRL AKITVDSLLK QGFMVKNIQK VTNEHPEDGF MLIPSYLFSG NDNLNLQIWI GHGRIFVHS
Source	Baculovirus
Target Names	At5g37930
Protein Names	Recommended name: E3 ubiquitin-protein ligase SINA-like 10 EC= 6.3.2.- Alternative name(s): Seven in absentia-like protein 10
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.