



# Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 56 (P0498E12.129, P0672D01.112)

<b>Product Code</b>	CSB-BP801037OFG
<b>Abbreviation</b>	Os08g0159800, LOC_Os08g06330
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q84UQ3
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDYTNAIHII PDAAGPDWA NAAAQGGDAG IWATEDDYNS QWNADGGGGG SSRAGSEQPP PGKKSRRGGG GEGGGNTSKS RAIGKMFFKT KLCCKFRAGT CPYVTNCNFA HGMEELRKPP PNWQEIVAAH EEATEAREEH QIPIMTSSGP TAGGDAGCGG GGGGGSGRAY KGRHCKKFYT DEGCPYGDAC TFLHDEQSKA RESVAISLSP SVGGGGGGGS YNSAAAAAAS ASAAAGNGPM QKPSNWKTRI CNKWEMTGYC PFGSKCHFAH GAAELHKYGG GLVDIDSRDA AATPDSKQAV VSAKAPAETA AASTTVLPHA DVYHLGVQAQ RSTIAGQRSQ QVQRPIQKWK GPKKISRIYG DWIDETE
<b>Source</b>	Baculovirus
<b>Target Names</b>	P0498E12.129, P0672D01.112
<b>Protein Names</b>	Recommended name: Zinc finger CCCH domain-containing protein 56 Short name= OsC3H56
<b>Expression Region</b>	1-367
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.