



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 56 (P0498E12.129, P0672D01.112)

Product Code	CSB-EP801037OFG-B
Abbreviation	Os08g0159800, LOC_Os08g06330
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84UQ3
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDYTNAIHII PDAAGPDWA NAAAQGGDAG IWATEDDYNS QWNADGGGGG SSRAGSEQPP PGKKSRRGGG GEGGGNTSKS RAIGKMFFKT KLCKFRAGT CPYVTNCNFA HGMEELRKP PNWQEIVAAH EEATEAREEH QIPIMTSSGP TAGGDAGCGG GGGGGSGRAY KGRHCKKFYT DEGCPYGDAC TFLHDEQSKA RESVAISLSP SVGGGGGGGS YNSAAAAAAS ASAAAGNGPM QKPSNWKTRI CNKWEMTGYC PFGSKCHFAH GAAELHKYGG GLVDIDSRDA AATPDSKQAV VSAKAPAETA AASTTVLPHA DVYHLGVQAQ RSTIAGQRSQ QVQRPIQKWK GPKKISRIYG DWIDETE
Source	E.coli
Target Names	P0498E12.129, P0672D01.112
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 56 Short name= OsC3H56
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.