



Recombinant Arabidopsis thaliana F-box protein At2g18780 (At2g18780)

Product Code	CSB-MP801038DOA
Abbreviation	At2g18780
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84V15
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQKQSTTMLS DLPKDVAEEV LSKFPVTSR RVRFTCKKLN TLSKDMSFTR KHIDKEEAKK KQSKEFHVVM MLDFRVSLFS LLLLNPPIE RLGQLVSLDG GAQVDISNIF HCEGFLCTT KDISRVVWN PCVGQTIWIK PRNSFNKWDR YALGYDVMKN HKVLRVDDG EYGRRHLLCE FEIYSFESDS WEVLVNPDW EIDFVHRGLS IDGNSYWFAR DKVVPYGGQD PSYFLLCDFD TTERFRPPLP LPFQPYFGDT VALSSVGDNQ IAVLCQRSSS PYTLKIWISS KVEPNAVSWN KLFLKVMKP LTGNSFSYSG GSFFVDEEKK MAVVLDKNIG GGFDPTRNIA YIIGKEGYFK KVDLGESRHL</p>
Source	Mammalian cell
Target Names	At2g18780
Protein Names	Recommended name: F-box protein At2g18780
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.