



# Recombinant Arabidopsis thaliana CBS domain-containing protein CBSX5 (CBSX5)

<b>Product Code</b>	CSB-EP801045DOA
<b>Abbreviation</b>	CBSX5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q84WQ5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MALSLLSYNV SDLCLGKPPPL RCLSSSSSSSV SDAIAALKSS EDTFLSVWNC NHDDDNTEC ECLGKISMAD VICHLSKDHD HSLCALNSSV SVLLPKTRSI VLHVQPSCSL IEAIDLIKQ AQNLIVPIHT KPYTKKKQHN DNVSVTTHH SNGQRFWIT QEDIIQFLLG FIAAFSPLPA MSLSDLGVIN STHTVVAVDY HSSASAVVSA VSNALAVQTS VAVVDGEGDD PFTSLIGEIS PMTLTCCDET AAAAVATLSA GDLMAYIDGA NPPELVQIV RNRLEDKGLI GLMSLFDSLS SYSTSSGYSS EEEAPVRTTS YGRSMSSSAR MARKSEAIVC NPKSSLMAVM IQAVAHRVNY AWWVEKDGCF VGMVTFVDIL KVFRKFLEND M
<b>Source</b>	E.coli
<b>Target Names</b>	CBSX5
<b>Protein Names</b>	Recommended name: CBS domain-containing protein CBSX5
<b>Expression Region</b>	1-391
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.