



Recombinant *Oryza sativa* subsp. japonica Zinc finger AN1 domain-containing stress-associated protein 14 (SAP14)

Product Code	CSB-MP801058OFG
Abbreviation	SAP14
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q852K8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MATKRKCPAN GDDGGVADLE PVAGGSFASP PPEKKAKLTV AVAVAVAPSS SSSATTAAAG EATAKREHGG FFAFARPENN TRLSVAVASS SSSASAAAEK AMAKLTVAGV APSSSASAAA AGKATAKREY GGFCAFARPD DKTRWRVAVA SSAAAAADAS YSSSSPATGE QPEANRCATC RRVKVLGTGFK CRCGGTFCGG HRYADEHGCG FDYKSSGREL IAKQNPVVVA DKLAFRI
Source	Mammalian cell
Target Names	SAP14
Protein Names	Recommended name: Zinc finger AN1 domain-containing stress-associated protein 14 Short name= OsSAP14
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.