



Recombinant Rice dwarf virus Minor outer capsid protein P9

Product Code	CSB-MP801061RCAO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q85440
Product Type	Recombinant Protein
Immunogen Species	Rice dwarf virus (isolate Fujian) (RDV)
Purity	>85% (SDS-PAGE)
Sequence	MGKIQDGLAI KRINDAITTF KNYKLGELKQ GGSMINTLS NVRAHVGLAW PAILRNCLIH TSSHLGFMKF MIDIATTWKV GAFTLLGSVG DEDPFTDVDL IYTKTCLHLG LKDNDLQFP EEFAYEANSF LEAQSMNAKV DMLTGVHNE DKYVFRMQSI SKFLKAYYTA SEDVAYLTGF IKPDDSKDSI LNAELLEAV TSEVLRVRNL ITTKIQKYIN LYEDSQLPHF RQAALSYIQD WVDVGGVPAA LPQPDTTDDE RPVTKPGPST PTVSKGVDEP EDEEMIRKKV ETSKDAPSKA DPPGNVSPRG VPALLEDDMS EMDMPDGFHD YLTREHENNF DLSQLGLAPS V
Source	Mammalian cell
Protein Names	Recommended name: Minor outer capsid protein P9
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.