



Recombinant Nascent polypeptide-associated complex subunit alpha (Y65B4BR.5)

Product Code	CSB-YP801217CXY
Abbreviation	Y65B4BR.5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86S66
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTGSTETRQK EVKEPQVDVS DDSNEAVEQ ELTEEQRRVA EAAGLGDHID KQAKQSRSEK KARKLFSKLG LKQVTGVS RV CIRKSKNILF VINKPDVFKS PGSDTYIIFG EAKIEDLTQH AQMSAIENLK PTREAPQLKT VEEDENEDVE FKEDSTGIEE KDIELVISQA NTRNKAIRA LKEADNDIVN AIMSLTM
Source	Yeast
Target Names	icd-2
Protein Names	Recommended name: Nascent polypeptide-associated complex subunit alpha Short name= NAC-alpha Alternative name(s): Alpha-NAC
Expression Region	1-197
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.