



Recombinant Human Serine/threonine-protein kinase tousled-like 2 (TLK2)

Product Code	CSB-YP801231HU
Abbreviation	TLK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UE8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MMEELHSLDP RRQELLEARF TGVGVSKGPL NSESSNQSLC SVGSLSDKEV ETPEKKQNDQ RNRKRKAEPY ETSQGKGTPR GHKISDYFEF AGGSAPGTSP GRSVPPVARS SPQHSLSNPL PRRVEQPLYG LDGSAAKEAT EEQSALPTLM SVMLAKPRLD TEQLAQRGAG LCFTFVSAQQ NSPSSTGSGN TEHSCSSQKQ ISIQHRQTQS DLTIEKISAL ENSKNSDLEK KEGRIDDLLR ANCDLRRQID EQQKMLEKYK ERLNRCVTMS KLLIEKSKQ EKMACRDKSM QDRLRLGHFT TVRHGASFTE QWTDGYAFQN LIKQQRERINS QREEIERQRK MLAKRKPPAM GQAPPATNEQ KQRKSKTNGA ENETPSSGNT ELKDTAPALG AHSLLRLTLA EYHEQEEIFK LRLGHLKKEE AEIQAELERL ERVRNLHIRE LKRIHNEDNS QFKDHPTLND RYLLHLHLLGR GGFSEVYKAF DLTEQRYVAV KIHQLNKNWR DEKKENYHKH ACREYRIHKE LDHPRIVKLY DYFSLDTSDF CTVLEYCEGN DLDFYLKQHK LMSEKEARSI IMQIVNALKY LNEIKPPIH YDLKPGNILL VNGTACGEIK ITDFGLSKIM DDDSYNSVDG MELTSQGAGT YWYLPPECFV VGKEPPKISN KVDVWSVGVI FYQCLYGRKP FGHNQSQQDI LQENTILKAT EVQFPKPVV TPEAKAFIRR CLAYRKEDRI DVQQLACDPY LLPHIRKSVS TSSPAGAAIA STSGASNSS SN
Source	Yeast
Target Names	TLK2
Protein Names	Recommended name: Serine/threonine-protein kinase tousled-like 2 EC=2.7.11.1 Alternative name(s): PKU-alpha Tousled-like kinase 2
Expression Region	1-772
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.