



Recombinant Human Tripartite motif-containing protein 74 (TRIM74)

Product Code	CSB-MP801235HU
Abbreviation	TRIM74
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UV6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAWQVSLLEL EDWLQCPICL EVFKESLMLQ CGHSYCKGCL VLSYHLDTK VRCPMCWQVV DGSSSLPNVS LAWVIEALRL PGDPEPKVCV HHRNPLSLFC EKDQELICGL CGLLGSHQHH PVTPVSTVCS RMKEELAALF SELKQEQKKV DELIAKLVKN RTRIVNESDV FSWVIRREFQ ELRHPVDEEK ARCLEGIGGH TRGLVASLDM QLEQAQGTRE RLAQAECVLE QFGNEDHHEF IWKFHSMASR
Source	Mammalian cell
Target Names	TRIM74
Protein Names	Recommended name: Tripartite motif-containing protein 74 Alternative name(s): Tripartite motif-containing protein 50C
Expression Region	1-250
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.